

NATIONAL RECALL ALERT CENTER - IMMEDIATE WARNING ALERT NOTIFICATION  
THE FOLLOWING HAVE BEEN RECALLED OR ARE SUBJECT TO FIELD CORRECTION  
THIS LISTING IS ISSUED BY - NATIONAL RECALL ALERT CENTER, WASHINGTON, D.C.

**THE FOLLOWING IS A CLASS I RECALL. THERE IS A REASONABLE PROBABILITY THAT THE USE OF THE AFFECTED PRODUCT MAY CAUSE SERIOUS INJURY OR DEATH IF THE PROBLEM IS NOT CORRECTED.**

**PRODUCT** Power Pro, Power ProMax, and MPower 1 Handpieces  
**RECALL NUMBER** Not available at the time of this publication  
**CODE** Catalog #s: PRO 5100, PRO 5100M, PRO 5200, PRO 5200M, PRO 5250, PRO 5250M, PRO 5300, PRO 5300M, PRO 5400, PRO 5400M, PRO 6200, PRO 6202, PRO6300, PRO 6400, PRO 6450  
 Please email [rhoda@recallalert.org](mailto:rhoda@recallalert.org) to have your serial numbers checked.  
**MANUFACTURED & RECALLED BY** Linvatec Corp, dba ConMed Linvatec, Largo, FL  
**QUANTITY** 14590 handpieces  
**DISTRIBUTION** Nationwide  
**REASON** The device has a switch problem where it may self-activate without pushing the trigger, continues running after releasing the trigger, and runs in unintended directions.  
**NOTE** Units manufactured and distributed after May 31, 2008 are not affected by this action.  
**CONTACT** ConMed Linvatec Customer Service at 800-237-0169

**THE FOLLOWING IS A CLASS I RECALL. THERE IS A REASONABLE PROBABILITY THAT THE USE OF THE AFFECTED PRODUCT MAY CAUSE SERIOUS INJURY OR DEATH IF THE PROBLEM IS NOT CORRECTED.**

**PRODUCT** Universal Cables, MC5057, manufactured from 1/3/2001 though 11/30/2006 and distributed from 1/24/2001 through 2/27/2009. The cable connects Linvatec MicroChoice Small Bone, MicroChoice Advantage and APEX Shaver System and Power Pro Electric I and II Large Bone System electric handpieces to the power source.  
**RECALL NUMBER** Not available at the time of this publication  
**CODE** Please email [rhoda@recallalert.org](mailto:rhoda@recallalert.org) to have your serial numbers checked.  
**MANUFACTURED & RECALLED BY** Linvatec Corp, dba ConMed Linvatec, Largo, FL  
**QUANTITY** 1768 cables  
**DISTRIBUTION** Nationwide  
**REASON** The cable may cause the electric handpiece to self-activate.  
**NOTE** Cables manufactured & distributed after December 1, 2006 at not affected by this action.  
**CONTACT** ConMed Linvatec Customer Service at 800-237-0169